

REAL-WORLD CASE STUDIES USING PURAPLY® AM

PuraPly AM® is an effective barrier to resist microbial colonization within the device and reduce microbes penetrating through the device.¹

Explore how PuraPly AM can take control of bioburden and biofilm re-formation to improve wound bed conditions and support healing in a wide range of wound types.¹⁻⁷



TABLE OF CONTENTS


TAKE CONTROL BEFORE BIOFILM RE-FORMS	2
BRIDGE BETWEEN DEBRIDEMENTS	3
TREATMENT ALGORITHM	4
OVERVIEW OF CASE STUDIES	5
CASE 1 Diabetic foot ulcer (DFU) with PuraPly® AM	6
CASE 2 Partial-thickness trauma wound with PuraPly AM	7
CASE 3 Heel pressure injury, which had failed a split-thickness skin graft, with PuraPly AM	8
CASE 4 Heel pressure injury with PuraPly AM	9
CASE 5 Surgical wound, which had failed a full-thickness skin graft, with PuraPly AM XT	10
CASE 6 Post-Mohs surgical wound (on right leg) with PuraPly AM XT transitioning to Affinity®	12
CASE 7 Post-Mohs surgical wound (on left leg) with PuraPly AM XT transitioning to Affinity	13
CASE 8 Venous leg ulcer (VLU) with PuraPly AM transitioning to Apligraf®	14
CASE 9 DFU with PuraPly AM transitioning to Apligraf	15
CASE 10 VLU with PuraPly AM transitioning to NuShield®	16
CASE 11 Trauma wound with PuraPly AM transitioning to NuShield	17

TAKE CONTROL BEFORE BIOFILM RE-FORMS

Use PuraPly® AM, an antimicrobial barrier, as early as day 1* to help move wounds out of the inflammatory phase.^{1,4,8,9}

**BIOFILM TRIGGERS NONSTOP INFLAMMATION
THAT STALLS WOUNDS⁹⁻¹¹**

**MANAGE BIOFILM TO HELP
CONTROL INFLAMMATION¹¹**



**BIOBURDEN &
BIOFILM**

**PROLONGED,
ELEVATED
INFLAMMATION**

EXCESS MMPs

**ECM
DEGRADATION**

Sharp debridement should be combined with an optimal barrier that contains a broad-spectrum antimicrobial and provides a sustained effect against bioburden and biofilm regrowth.¹¹

*Certain local coverage determinations may not allow skin substitute use until after conventional care at week 4.
MMPs = matrix metalloproteinases; ECM = extracellular matrix

BRIDGE BETWEEN DEBRIDEMENTS

Only PuraPly® AM, an antimicrobial barrier, combines native, cross-linked ECM plus broad-spectrum PHMB.¹



Scan to learn more about
how PuraPly AM works.

**Provides a sustained
antimicrobial barrier effect³**

**Eliminates the need for at-home
primary dressing changes^{1,2}**




**Helps prevent biofilm
re-formation³**

**Keeps you in control of
the healing environment**

NATIVE, CROSS-LINKED ECM

Enhances resistance
to enzymatic
degradation^{12,13}




BROAD-SPECTRUM PHMB

Helps control
bioburden³⁻⁶

TREATMENT ALGORITHM

Organogenesis—Elevating the standard of care with a broad portfolio of innovative solutions



* E = **Exudate** well-controlled; I = No signs of clinical **Infection**; G = Healthy **Granulation** tissue present

† Please refer to the Apligraf Package Insert for complete prescribing information and contraindications.

VLU = venous leg ulcer; DFU = diabetic foot ulcer

OVERVIEW OF CASE STUDIES

Patient demographics and wound specifics


PATIENTS (n=10)	
MALES	4
FEMALES	6
MEAN AGE (YEARS)	73
AGE RANGE (YEARS)	53-96

BASELINE WOUND AREA	
MEAN	20.7 cm ²
MEDIAN	12.0 cm ²
RANGE	3.8-80.5 cm ²

BASELINE WOUND DURATION



Types of wounds (n=11)



- 2** DIABETIC FOOT ULCERS (DFUs)
- 2** POST-MOHS SURGICAL WOUNDS*
- 2** PRESSURE INJURIES
- 1** SURGICAL WOUND
- 2** TRAUMA WOUNDS
- 2** VENOUS LEG ULCERS (VLUs)

*One patient had two post-Mohs surgical wounds.

PuraPly® AM, an antimicrobial barrier, supported healing of a DFU at 6 weeks following 5 applications

71-year-old female with a DFU resulting from a 2nd degree burn (due to boiling water) on the left heel

WOUND SPECIFICS

- Wound area was 14.8 cm² at the time of PuraPly AM initiation
- Wound was present for 13 weeks (3 months)

PATIENT HISTORY

- Wound healing was complicated by diabetes, excessive pressure, and peripheral neuropathy
- Chronic kidney disease, peripheral vascular disease, hypertension, and severe rheumatoid arthritis
- Surgical history included OR debridement and irrigating NPWT; patient was previously hospitalized

PREVIOUS TREATMENTS

Enzymatic debridement, NPWT, and Dermagraft® (1 application)

PURAPLY AM APPLICATION SCHEDULE

- Weekly applications for 5 weeks
- Offloading boot was used throughout the course of treatment

1st PuraPly AM application



WEEK 1

Pre-debridement

Date: 12/16/15

Wound Size: 3.8 x 3.9 x 0.5 cm

Wound Area: 14.8 cm²

2nd PuraPly AM application



WEEK 2

Pre-debridement

Date: 12/23/15

Wound Size: 2.5 x 3.3 x 0.1 cm

Wound Area: 8.3 cm²

3rd PuraPly AM application



WEEK 3

Pre-debridement

Date: 12/30/15

Wound Size: 1.8 x 2.9 x 0.1 cm

Wound Area: 5.2 cm²

4th PuraPly AM application



WEEK 4

Pre-debridement

Date: 1/6/16

Wound Size: 1.3 x 2.0 x 0.1 cm

Wound Area: 2.6 cm²

5th PuraPly AM application



WEEK 5

Pre-debridement

Date: 1/13/16

Wound Size: 0.5 x 1.0 x 0.1 cm

Wound Area: 0.5 cm²

Complete wound closure



WEEK 6

Date: 1/20/16

Wound remained closed



WEEK 12

Date: 2/29/16

PuraPly® AM, an antimicrobial barrier, supported healing of a partial-thickness trauma wound at 4 weeks following one application

53-year-old female with a trauma wound resulting from a hematoma on the left shin

WOUND SPECIFICS

- Wound area was 80.5 cm² at the time of PuraPly AM initiation
- Partial-thickness wound with islands of epithelium noted
- Wound was present for 4 weeks

PATIENT HISTORY

- Wound healing was complicated by venous insufficiency
- COPD, peripheral vascular disease, hypertension, diabetes, and obesity
- Surgical history included bilateral vein ablation

PREVIOUS TREATMENTS

Four weeks of compression and enzymatic debridement

PURAPLY AM APPLICATION SCHEDULE

- 1 application
- Compression with Kerlix and Coban was used throughout course of treatment

Initial presentation



4 weeks prior to PuraPly AM

Post-debridement

Date: 6/9/16

Wound Size: 12.0 x 7.5 x 0.0 cm

Wound Area: 90.0 cm²

No PuraPly AM application



1 week prior to PuraPly AM

Post-debridement

Date: 6/23/16

Wound Size: 12.0 x 7.0 x 0.0 cm

Wound Area: 84.0 cm²

1st PuraPly AM application



WEEK 1

Pre-debridement

Date: 6/30/16

Wound Size: 11.5 x 7.0 x 0.0 cm

Wound Area: 80.5 cm²

No PuraPly AM application*



WEEK 2

Post-debridement

Date: 7/7/16

Wound Size: 0.8 x 0.4 x 0.0 cm

Wound Area: 0.3 cm²

No PuraPly AM application†



WEEK 3

Post-debridement

Date: 7/14/16

Wound Size: 1.5 x 2.5 x 0.0 cm

Wound Area: 3.8 cm²

Complete wound closure



WEEK 4

Date: 7/20/16

*Dramatic reduction in wound size with only minimal open areas throughout the wound bed. Physician decided PuraPly AM application was not needed.

†Devitalized tissue removed revealing a healthy wound bed. Physician decided PuraPly AM application was not needed.

PuraPly® AM, an antimicrobial barrier, supported healing of a heel pressure injury at 10 weeks following 9 applications, after failing a split-thickness skin graft (STSG)

59-year-old male with a pressure injury on the left heel

WOUND SPECIFICS

- Wound area was 12.0 cm² at the time of PuraPly AM initiation
- STSG failed 11 days prior to PuraPly AM application

PATIENT HISTORY

- Left heel pressure injury developed during hospitalization
- Diabetes, peripheral arterial disease, osteomyelitis, and edema
- Surgical history included right BKA

PREVIOUS TREATMENTS

Debridement, partial calcaneotomy, tissue expander removal, NPWT, and failed STSG

PURAPLY AM APPLICATION SCHEDULE

- Weekly applications for 9 weeks
- Offloading with a customized cut out was used throughout the course of treatment


1st PuraPly AM application



WEEK 1

Pre-debridement showing failed STSG
Date: 6/29/15
Wound Size: 2.8 x 4.3 x 0.1 cm
Wound Area: 12.0 cm²

2nd PuraPly AM application



WEEK 2

Pre-debridement
Date: 7/6/15
Wound Size: 2.6 x 4.1 x 0.1 cm
Wound Area: 10.7 cm²


3rd PuraPly AM application



WEEK 3

Pre-debridement
Date: 7/13/15
Wound Size: 2.0 x 3.6 x 0.1 cm
Wound Area: 7.2 cm²


4th PuraPly AM application



WEEK 4

Pre-debridement
Date: 7/20/15
Wound Size: 2.0 x 3.3 x 0.1 cm
Wound Area: 6.6 cm²


5th PuraPly AM application



WEEK 5

Pre-debridement
Date: 7/27/15
Wound Size: 1.4 x 1.8 x 0.1 cm
Wound Area: 2.5 cm²


6th PuraPly AM application



WEEK 6

Pre-debridement
Date: 8/3/15
Wound Size: 1.2 x 1.5 x 0.1 cm
Wound Area: 1.8 cm²


7th PuraPly AM application



WEEK 7

Pre-debridement
Date: 8/10/15
Wound Size: 0.6 x 0.6 x 0.1 cm
Wound Area: 0.4 cm²


8th PuraPly AM application



WEEK 8

Pre-debridement
Date: 8/17/15
Wound Size: 0.6 x 0.6 x 0.1 cm
Wound Area: 0.4 cm²

9th PuraPly AM application



WEEK 9

Pre-debridement
Date: 8/24/15
Wound Size: 0.4 x 0.3 x 0.1 cm
Wound Area: 0.1 cm²

Complete wound closure



WEEK 10

Date: 8/31/15
Patient ready to be molded for crow boot and custom AFO.

PuraPly® AM, an antimicrobial barrier, supported healing of a heel pressure injury at 10 weeks following 9 applications

67-year-old male with a pressure injury on the left heel at the site of a previously closed wound

WOUND SPECIFICS

- Wound area was 18.0 cm² at the time of PuraPly AM initiation
- Wound recurrence was likely due to crow boot shearing while walking on right prosthesis
- Wound was present for 9 weeks

PATIENT HISTORY

- Wound healing complicated by diabetes, neuropathy, and peripheral vascular disease
- Hypertension, gout, ESRD, hyperlipidemia, anemia, and osteomyelitis (right heel)
- Surgical history included partial calcaneotomy bilaterally, right BKA, and surgical resection of left heel

PREVIOUS TREATMENTS

NPWT

PURAPLY AM APPLICATION SCHEDULE

- Weekly applications for 9 weeks
- Offloading with total contact cast was used throughout the course of treatment

1st PuraPly AM application




WEEK 1

Pre-debridement

Date: 6/29/15

Wound Size: 4.0 x 4.5 x 0.2 cm
Wound Area: 18.0 cm²

2nd PuraPly AM application




WEEK 2

Pre-debridement

Date: 7/6/15

Wound Size: 3.7 x 4.2 x 0.2 cm
Wound Area: 15.5 cm²

3rd PuraPly AM application




WEEK 3

Pre-debridement

Date: 7/13/15

Wound Size: 2.4 x 4.0 x 0.2 cm
Wound Area: 9.6 cm²

4th PuraPly AM application




WEEK 4

Pre-debridement

Date: 7/20/15

Wound Size: 2.3 x 4.0 x 0.2 cm
Wound Area: 9.3 cm²

5th PuraPly AM application




WEEK 5

Pre-debridement

Date: 7/27/15

Wound Size: 2.1 x 3.3 x 0.1 cm
Wound Area: 6.9 cm²

6th PuraPly AM application




WEEK 6

Pre-debridement

Date: 8/3/15

Wound Size: 1.4 x 3.0 x 0.1 cm
Wound Area: 4.2 cm²

7th PuraPly AM application




WEEK 7

Pre-debridement

Date: 8/10/15

Wound Size: 1.2 x 2.5 x 0.1 cm
Wound Area: 3.0 cm²

8th PuraPly AM application



WEEK 8

Pre-debridement

Date: 8/17/15

Wound Size: 0.7 x 2.1 x 0.1 cm
Wound Area: 1.5 cm²

9th PuraPly AM application




WEEK 9

Pre-debridement

Date: 8/24/15

Wound Size: 0.2 x 0.7 x 0.1 cm
Wound Area: 0.1 cm²

Complete wound closure



WEEK 10

Date: 8/31/15

PuraPly® AM XT, an antimicrobial barrier, supported healing of a surgical wound at 6 weeks following 5 applications, after failing a full-thickness skin graft (FTSG)

79-year-old female with a surgical wound on the left medial leg

WOUND SPECIFICS

- Wound area was 12.0 cm² at the time of PuraPly AM XT initiation
- Wound was present for 6 weeks

PATIENT HISTORY

- Melanoma and cardiac ablation
- Surgical history included melanoma excision, which was closed using a FTSG

PREVIOUS TREATMENTS

Xeroform and Unna's boot

PURAPLY AM XT APPLICATION SCHEDULE

- Weekly applications for 5 weeks

Initial presentation



1 week prior to PuraPly AM XT

Post-debridement
Date: 12/2/20
Wound Size: 3.1 x 4.3 x 0.3 cm
Wound Area: 13.3 cm²

1st PuraPly AM XT application



WEEK 1

Pre-debridement
Date: 12/9/20
Wound Size: 3.0 x 4.0 x 0.3 cm
Wound Area: 12.0 cm²

2nd PuraPly AM XT application



WEEK 2

Post-debridement
Date: 12/16/20
Wound Size: 2.5 x 3.9 x 0.3 cm
Wound Area: 9.8 cm²

3rd PuraPly AM XT application



WEEK 3

Post-debridement
Date: 12/23/20
Wound Size: 2.0 x 3.0 x 0.3 cm
Wound Area: 6.0 cm²

4th PuraPly AM XT application



WEEK 4

Pre-debridement
Date: 12/30/20
Wound Size: 2.0 x 3.0 x 0.3 cm
Wound Area: 6.0 cm²

5th PuraPly AM XT application



WEEK 5

Post-debridement
Date: 1/6/21
Wound Size: 1.3 x 0.9 x 0.3 cm
Wound Area: 1.2 cm²

Complete wound closure



WEEK 6

Date: 1/13/21

Thin layer of epithelium present.

**CASE STUDIES
UTILIZING PURAPLY® AM
AND TRANSITIONING
TO ANOTHER
ORGANOGENESIS SOLUTION**

PuraPly® AM XT, an antimicrobial barrier, transitioning to Affinity®, a fresh amniotic membrane wound covering, supported healing of a post-Mohs surgical wound at 8 weeks

95-year-old female with a surgical wound on the right leg (same patient as on page 13)

WOUND SPECIFICS

- Wound area was 10.5 cm² at the time of PuraPly AM XT initiation
- Wound was present for 1 week

PATIENT HISTORY

- Heart murmur, total hip replacement, and venous insufficiency
- Surgical history included Mohs surgery to remove a squamous cell carcinoma on the right leg

PREVIOUS TREATMENTS

Xeroform and 2-layer compression

PURAPLY AM XT APPLICATION SCHEDULE

- Weekly applications for 6 weeks

Initial presentation



1 week prior to PuraPly AM XT

Post-debridement
Date: 5/26/21
Wound Size: 2.8 x 2.8 x 0.4 cm
Wound Area: 7.8 cm²

1st PuraPly AM XT application



WEEK 1

Post-debridement
Date: 6/2/21
Wound Size: 3.0 x 3.5 x 0.5 cm
Wound Area: 10.5 cm²


2nd PuraPly AM XT application



WEEK 2

Post-debridement
Date: 6/11/21
Wound Size: 2.2 x 2.5 x 0.5 cm
Wound Area: 5.5 cm²


3rd PuraPly AM XT application



WEEK 3

Post-debridement
Date: 6/16/21
Wound Size: 2.2 x 2.5 x 0.5 cm
Wound Area: 5.5 cm²


4th PuraPly AM XT application



WEEK 4

Post-debridement
Date: 6/23/21
Wound Size: 1.7 x 2.3 x 0.5 cm
Wound Area: 3.9 cm²

5th PuraPly AM XT application



WEEK 5

Post-debridement
Date: 6/30/21
Wound Size: 1.2 x 1.4 x 0.5 cm
Wound Area: 1.7 cm²

6th PuraPly AM XT application



WEEK 6

Post-debridement
Date: 7/7/21
Wound Size: 1.2 x 1.4 x 0.5 cm
Wound Area: 1.7 cm²

Transition to Affinity*



WEEK 7

Post-debridement
Date: 7/14/21
Wound Size: 0.4 x 0.2 x 0.2 cm
Wound Area: 0.08 cm²

Complete wound closure



WEEK 8

7/21/21
Patient received one application of Affinity.

*Reduction in wound size with healthy granulation tissue and less slough. Clinician transitioned to Affinity.

PuraPly® AM XT, an antimicrobial barrier, transitioning to Affinity®, a fresh amniotic membrane wound covering, supported healing of a post-Mohs surgical wound at 7 weeks

95-year-old female with a surgical wound on the left leg (same patient as on page 12)

WOUND SPECIFICS

- Wound area was 3.8 cm² at the time of PuraPly AM XT initiation
- Wound was present for 4 weeks

PATIENT HISTORY

- Heart murmur, total hip replacement, and venous insufficiency
- Surgical history included Mohs surgery to remove a basal cell carcinoma on the left leg
- Patient underwent this procedure while she was being treated for another surgical wound on the right leg

PREVIOUS TREATMENTS

Xeroform and oral antibiotic treatment prior to PuraPly AM XT application

PURAPLY AM XT APPLICATION SCHEDULE

- Weekly applications for 2 weeks

Initial presentation



3 weeks prior to PuraPly AM XT

Post-debridement
Date: 6/23/21
Wound Size: 2.8 x 2.7 x 0.3 cm
Wound Area: 7.6 cm²

1st PuraPly AM XT application



WEEK 1

Post-debridement
Date: 7/14/21
Wound Size: 1.9 x 2.0 x 0.7 cm
Wound Area: 3.8 cm²


2nd PuraPly AM XT application



WEEK 2

Post-debridement
Date: 7/21/21
Wound Size: 1.4 x 1.7 x 0.6 cm
Wound Area: 2.4 cm²


Transition to Affinity*



WEEK 3

Post-debridement
Date: 7/28/21
Wound Size: 1.1 x 1.1 x 0.6 cm
Wound Area: 1.2 cm²


2nd Affinity application



WEEK 4

Post-debridement
Date: 8/4/21
Wound Size: 0.6 x 0.9 x 0.4 cm
Wound Area: 0.5 cm²

3rd Affinity application



WEEK 5

Post-debridement
Date: 8/11/21
Wound Size: 0.5 x 0.8 x 0.4 cm
Wound Area: 0.4 cm²

No product application required



WEEK 6

Post-debridement
Date: 8/18/21
Wound Size: 0.1 x 0.1 x 0.1 cm
Wound Area: 0.01 cm²

Complete wound closure



WEEK 7

Date: 8/25/21

*Reduction in wound size with healthy granulation tissue and less slough. Clinician transitioned to Affinity.

PuraPly® AM, an antimicrobial barrier, transitioning to Apligraf®, bioengineered with living cells, supported healing of a VLU at 12 weeks

58-year-old female with a VLU on the distal lateral aspect of the left leg

WOUND SPECIFICS

- Wound area was 35.0 cm² at the time of PuraPly AM initiation
- Wound was present for approximately 13 weeks (3+ months)

PATIENT HISTORY

- Congestive heart failure, peripheral vascular disease, hypertension, and diabetes

PREVIOUS TREATMENTS

Wet-to-dry dressings and Unna's boot compression

PURAPLY AM APPLICATION SCHEDULE

- Weekly applications for 5 weeks
- Compression with Kerlix and Coban used during PuraPly AM applications

Initial presentation



1 week prior to PuraPly AM

No debridement*


Date: 5/12/16

Wound Size: 6.0 x 4.5 x 0.5 cm

Wound Area: 27.0 cm²

*Patient could not tolerate debridement due to pain.
OR debridement scheduled for the following week.

1st PuraPly AM application



WEEK 1

Post-debridement

Date: 5/18/16

Wound Size: 7.0 x 5.0 x 0.5 cm

Wound Area: 35.0 cm²

2nd PuraPly AM application



WEEK 2


Post-debridement

Date: 5/26/16

Wound Size: 6.5 x 5.0 x 0.5 cm

Wound Area: 32.5 cm²

3rd PuraPly AM application



WEEK 3

Post-debridement

Date: 6/3/16

Wound Size: 6.0 x 5.0 x 0.5 cm

Wound Area: 30.0 cm²

4th PuraPly AM application



WEEK 4

Post-debridement

Date: 6/7/16

Wound Size: 5.0 x 4.5 x 0.5 cm

Wound Area: 22.5 cm²

Transition to Apligraf†



WEEK 5

Pre-debridement

Date: 6/16/16

Wound Size: 5.0 x 3.5 x 0.5 cm

Wound Area: 17.5 cm²

Complete wound closure



WEEK 12

Date: 8/2/16

Wound closed following 3 applications of Apligraf.

†Bioburden well-managed and reduction in hypergranulation, malodour, and drainage. Clinician transitioned to Apligraf to stimulate healing.

PuraPly® AM, an antimicrobial barrier, transitioning to Apligraf®, bioengineered with living cells, supported healing of a DFU at 13 weeks

80-year-old male with a Wagner Grade 2 DFU on the medial aspect of the 1st left metatarsal head

WOUND SPECIFICS

- Wound area was 10.0 cm² at the time of PuraPly AM initiation
- Wound was present 44 weeks (9+ months)

PATIENT HISTORY

- Wound healing was complicated by diabetes and peripheral vascular disease
- Atrial fibrillation, chronic obstructive pulmonary disease, glaucoma, hypertension, and iron deficiency anemia

PREVIOUS TREATMENTS

Weekly debridements and various dressings

PURAPLY AM APPLICATION SCHEDULE

- Weekly applications for 5 weeks
- Offloading shoe was used throughout the course of treatment

1st PuraPly AM application



WEEK 1

Post-debridement

Date: 9/21/15

Wound Size: 4.0 x 2.5 x 0.2 cm

Wound Area: 10.0 cm²

2nd PuraPly AM application



WEEK 2

Post-debridement

Date: 9/28/15

Wound Size: 3.8 x 2.1 x 0.2 cm

Wound Area: 8.0 cm²

3rd PuraPly AM application



WEEK 3

Post-application

Date: 10/5/15

Wound Size: 3.6 x 2.5 x 0.2 cm

Wound Area: 9.0 cm²

4th PuraPly AM application



WEEK 4

Post-debridement

Date: 10/12/15

Wound Size: 3.6 x 1.3 x 0.2 cm

Wound Area: 4.7 cm²

5th PuraPly AM application



WEEK 5

Post-debridement

Date: 10/19/15

Wound Size: 3.6 x 1.9 x 0.2 cm

Wound Area: 6.8 cm²

Transition to Apligraf*



WEEK 6

Post-debridement

Date: 10/26/15

Wound Size: 2.6 x 1.3 x 0.1 cm

Wound Area: 3.4 cm²

Complete wound closure



WEEK 13

Date: 12/16/15

Wound closed following 5 applications of Apligraf.

*Reduction in wound size with bioburden well-managed and healthy granulation tissue present. Clinician transitioned to Apligraf to stimulate healing.

PuraPly® AM, an antimicrobial barrier, transitioning to NuShield®, a complete dehydrated placental allograft wound covering, supported healing of a VLU at 11 weeks

96-year-old male with a VLU on the left anterior shin

WOUND SPECIFICS

- Wound area was 5.7 cm² at the time of PuraPly AM initiation
- Wound was present for 17 weeks (4 months)

PATIENT HISTORY

- Diabetes, hypertension, neuropathy, and venous insufficiency

PREVIOUS TREATMENTS

Enzymatic debridement, Medihoney®, and Promogran Prisma®

PURAPLY AM APPLICATION SCHEDULE

- Weekly applications for 3 weeks
- Compression was not utilized; patient was asked to elevate

1st PuraPly AM application



Post-debridement
Wound Size: 2.7 x 2.1 x 0.1 cm
Wound Area: 5.7 cm²

2nd PuraPly AM application




Pre-debridement
Wound Size: 2.7 x 2.1 x 0.1 cm
Wound Area: 5.7 cm²

3rd PuraPly AM application



Pre-debridement
Wound Size: 2.5 x 1.9 x 0.1 cm
Wound Area: 4.8 cm²


Transition to NuShield*



Post-debridement
Wound Size: 2.1 x 1.5 x 0.1 cm
Wound Area: 3.2 cm²

*Bioburden well-managed with an increase in granulation tissue. Clinician transitioned to NuShield.

Complete wound closure



Wound closed following 4 applications of NuShield.

Wound remained closed



PuraPly® AM, an antimicrobial barrier, transitioning to NuShield®, a complete dehydrated placental allograft wound covering, supported healing of a trauma wound at 16 weeks

69-year-old female with a trauma wound on the right lateral leg complicated by immunosuppressive therapy

WOUND SPECIFICS

- Wound was initially caused by trauma, was sutured but dehiscd and became infected
- Wound area was 25.5 cm² at the time of PuraPly AM initiation
- Wound was present for 4 weeks

PATIENT HISTORY

- Wound healing was complicated by immunosuppressive therapy and venous insufficiency
- Patient was on several medications to treat rheumatoid arthritis including etanercept, leflunomide, and prednisone
- Hypertension, hyperthyroidism, and GERD

PREVIOUS TREATMENTS

Enzymatic debridement; oral antibiotics were used to treat the wound infection prior to PuraPly AM initiation

PURAPLY AM APPLICATION SCHEDULE

- Weekly applications for 6 weeks
- Compression with Tubigrip was used during PuraPly AM applications

1st PuraPly AM application



WEEK 1
Pre-debridement
Date: 9/3/20
Wound Size: 8.5 x 3.0 x 0.5 cm
Wound Area: 25.5 cm²

2nd PuraPly AM application



WEEK 2
Pre-debridement
Date: 9/11/20
Wound Size: 8.3 x 3.4 x 0.5 cm
Wound Area: 28.2 cm²

3rd PuraPly AM application



WEEK 3
Pre-debridement
Date: 9/16/20
Wound Size: 8.5 x 3.2 x 0.4 cm
Wound Area: 27.2 cm²

4th PuraPly AM application



WEEK 4
Pre-debridement
Date: 9/25/20
Wound Size: 7.8 x 3.0 x 0.4 cm
Wound Area: 23.4 cm²

5th PuraPly AM application



WEEK 5
Pre-debridement
Date: 10/5/20
Wound Size: 7.0 x 3.0 x 0.1 cm
Wound Area: 21.0 cm²

6th PuraPly AM application



WEEK 6
Pre-debridement
Date: 10/12/20
Wound Size: 6.8 x 2.6 x 0.1 cm
Wound Area: 17.7 cm²

Transition to NuShield®



WEEK 7
Pre-debridement
Date: 10/19/20
Wound Size: 6.2 x 2.5 x 0.1 cm
Wound Area: 15.5 cm²

Complete wound closure



WEEK 16
Date: 12/21/20
Wound closed following 5 applications of NuShield.

*Bioburden was well-managed. Clinician transitioned to NuShield.

PURAPLY® AM KEEPS YOU IN CONTROL



Give your patients the new standard of care.

Control bioburden right from the start and all week long with PuraPly AM, an effective barrier to resist microbial colonization within the device and reduce microbes penetrating through the device.¹⁻⁷

Bioburden control as early as day 1 and all week long

Following sharp debridement, applying PuraPly AM as an antimicrobial barrier provides a foundation of care for weeklong control of bioburden and prevention of biofilm re-formation.¹⁻⁷

The unique power of plus for sustained control

Native, cross-linked ECM plus broad-spectrum PHMB acts as a bridge between weekly debridements by producing a sustained antimicrobial barrier effect.¹⁻⁴

Expanding evidence base of real-world clinical and scientific data

PuraPly AM gives you the confidence of control, with expanding clinical evidence in a variety of acute and chronic wounds and scientific data.^{3,7,18}

References: **1.** PuraPly Antimicrobial [package insert]. Canton, MA: Organogenesis Inc; 2020. **2.** PuraPly XT Extra Fenestrated [package insert]. Canton, MA: Organogenesis Inc. 2020. **3.** Davis SC, et al. *Int Wound J*. 2021;1-14. <https://doi.org/10.1111/iwj.13600>. **4.** Brantley J, et al. *Wounds Int*. 2016;7(3):1-5. **5.** Gilbert P, et al. *J Appl Microbiol*. 2005;99(4):703-715. **6.** Hübner NO, et al. *Skin Pharmacol Physiol*. 2010;23(1 suppl):17-27. **7.** Bain MA, et al. *J Comp Eff Res*. 2020;9(10):691-703. **8.** Carpenter S, et al. *Wounds*. 2016;28(6 suppl):S1-S20. **9.** Frykberg RG, Banks J. *Adv Wound Care*. 2015;4(9):560-582. **10.** Gibson D, et al. *Wounds Int*. 2009;1(1):1-6. **11.** Schultz G, et al. *Wound Repair Regen*. 2017;25(5):744-757. **12.** Khor E. *Biomaterials*. 1997;18(2):95-105. **13.** Billiar K, et al. *J Biomed Mater Res*. 2001;56(1):101-108. **14.** Phillips TJ, et al. *J Am Acad Dermatol*. 2000;43(4):627-630. **15.** Gelfand JM, et al. *J Invest Dermatol*. 2002;119(6):1420-1425. **16.** Sheehan P, et al. *Diabetes Care*. 2003;26(6):1879-1882. **17.** Attinger CE, et al. *Plast Reconstr Surg*. 2006;117(Suppl):72S-109S. **18.** Oropallo AR. *Plast Reconstr Surg Glob Open*. 2019;7:e2047.